



Analysis of Radioactive Materials

Capitalizing on 50+ years of success in developing and servicing shielded microanalytical instruments for nuclear fuel characterization, irradiated materials behavior investigation and radioactive waste management, CAMECA now launches the **SKAPHIA Shielded Electron Probe Initiative**, bringing together key players in nuclear research & industry towards the development of the next generation Shielded EPMA.

A long success story in supporting nuclear science & industry


Since the installation of the first shielded Electron Probe Microanalyzer in 1978, CAMECA has developed customized microprobes for major post-irradiation examination facilities such as Idaho National Lab (USA), Korea Atomic Energy Research Institute, SCK-CEN (Belgium), Institute for Transuranium Elements (Germany), EDF and several CEA research centers in France (cf. iaea.org).

Over the years, we have perfected our EPMA instrumentation in order to deliver ever higher spatial resolution and better detection limits, while also developing shielding methods and systems allowing nuclear scientists to efficiently analyze radioactive materials.

The SKAPHIA Shielded Electron Probe Initiative

Today, after five generations of CAMECA shielded electron microprobe analyzers, among which several SX50-R are still in operation, CAMECA is ready to develop its next generation Shielded EPMA.


While providing a safe environment for manipulating and analyzing nuclear samples, SKAPHIA will provide benchmark analytical performance. It will allow scientists to gain a deeper understanding of fuel performance, to explore irradiated material behavior and radiation damage processes, to develop innovative alloys and structural materials, to optimize the nuclear fuel cycle and to achieve better nuclear waste management, thus contributing to a safer world.



SKAPHIA will offer benchmark EPMA analytical performance.

Derived from CAMECA state-of-the-art EPMA instruments, our next generation Shielded EPMA will analyze almost all elements of the periodic table, revealing compositional information for both major and trace elements of radioactive sample. The information will be obtained from sub-micron areas with ultimate precision and accuracy.

The system will be designed to accept samples emitting β and γ radiations with a maximum acceptable γ radiation level of 111 GBq at an energy of 0.75 MeV.




- Column, spectrometers and sample stage will be installed in a hot cell (lead or concrete shielded room).
- Remote-manipulators and/or ball manipulators will be used to insert and mount the radioactive samples.
- The instrument will be fully remote-controlled (stage, column, diaphragms...) with electronics and computer deported outside the "hot" cell environment.
- All WDS analyzers and detectors will be shielded to prevent the background caused by the γ radiations.
- The stage will be made of Denal material.
- The Secondary Electron detector will have a special orientation to avoid γ ray interference.
- Dedicated CAMECA EPMA software will provide all necessary features for quantitative analysis, X ray mapping, line profile acquisition and data processing.

Applications

CAMECA Shielded EPMA have been used in a wide range of nuclear applications: study of the radiation damage processes in materials and complex micro-structures, analysis of materials for the nuclear fuel cycle, R&D of advanced alloys and innovative structural materials, nuclear reactor maintenance optimization, nuclear waste management...

Shown here are results from EPMA analyses of MgO pellet including AmOx particules irradiated up to 318 effective full power days (EFPD) in the Phénix reactor. The results demonstrated the Am transmutation rate reached was close to 95.5%.

Courtesy of Karine Hanifi and Jérôme Lamontagne, CEA Cadarache, France.





SKAPHIA

SHIELDED ELECTRON PROBE INITIATIVE

Analysis of radioactive materials
Post-irradiation examinations of nuclear fuel
Characterization of advanced alloys
Nuclear waste management



Selected bibliography

Electron Probe Microanalysis of Transuranium Elements in Nuclear Fuels and Waste Matrixes. M. Tribet, X. Deschanel, S. Peugeot, C. Jegou, K. Hanifi, J. Noirot, J. Lamontagne. Presentation at EMAS Conference, Marcoule, France, May 2016.

Magay Heat Capacity of Bi₂UO₆. K. Popa, O. Beneš, P. E. Raison, J-C. Griveau, P. Pöml, E. Colineau, R.J.M. Konings, J. Somers. *Journal of Nuclear Materials*, Vol. 465, p. 653-656, 2015.

ECRIX-H Irradiation: Post-Irradiation Examinations and Simulations. S. Béjaoui, J. Lamontagne, E. Esbelin, J.M. Bonnerot, E. Brunon, P. Bourdot, Y. Pontillon. Presentation at FP7 FAIRFUELS Workshop, Stockholm, Sweden, February 2011.

Chemical States of Fission Products and Actinides in Irradiated Oxide Fuels Analyzed by Thermodynamic Calculation and Post-Irradiation Examination. K. Kurosaki, K. Tanaka, M. Osaka, Y. Ohishi, H. Muta, M. Uno, S. Yamanaka. *Progress in Nuclear Science and Technology*, Vol. 2, p.5-8, 2011.

Microstructural Evolution and Am Migration Behavior in Am-Containing MOX Fuels at the Initial Stage of Irradiation. K. Tanaka, S. Miwa, I. Sato, M. Osaka, T. Hirose, H. Obayashi, S. Koyama, H. Yoshimochi, K. Tanaka. Presentation at the 10th OECD Nuclear Energy Agency Information Exchange Meeting on Actinide and Fission Product Partitioning and Transmutation, Mito, Japan, October 2008.

On the Oxidation State of UO₂ Nuclear Fuel at a Burn-Up of Around 100 MWd/kgHM. C.T. Walker, V.V. Rondinella, D. Papaioannou, S. Van Winckel, W. Goll, R. Manzel. *Journal of Nuclear Materials*, Vol. 345, p. 192-205, 2005.

Analysis of High Radioactive Materials in Irradiated DUPIC SIMFUEL Using EPMA. Jung, Yang Hong, Yoo, Bang Ok, Joo Yong Sun, Kim, Hee Mun, Jung In Ha, Kim, Myung Han. *Journal of the Korean Radioactive Waste Society*, Vol. 2(2), p. 125-133, 2004.

Multiple Voltage Electron Probe Microanalysis of Fission Gas Bubbles in Irradiated Nuclear Fuel. M. Verwerft. *Journal of Nuclear Materials*, Vol. 282, p. 97-111, 2000.

CAMECA is the world premier provider of microanalytical instrumentation. We deliver cutting-edge science and metrology solutions, and offer our customers unparalleled support and maintenance service through the comprehensive AMECARE program.

Corporate Headquarters:

29 Quai des Grésillons • 92622 Gennevilliers Cedex • FRANCE • Phone: +33 1 43 34 62 00 • cameca.info@ametek.com

CAMECA Subsidiaries or International AMETEK Offices:

AMETEK do Brasil CAMECA Business Unit • cameca@ametek.com.br

AMETEK China CAMECA Business Unit • cameca-china.sales@ametek.com.cn

CAMECA Germany • sales.cameca@ametek.de

AMETEK India CAMECA Business Unit • cameca-india@ametek.com

AMETEK Japan CAMECA Business Unit • cameca-japan.info@ametek.co.jp

AMETEK Korea CAMECA Business Unit • cameca-korea.info@ametek.co.kr

CAMECA Taiwan Corp. Ltd. • sales.cameca@ametek.com.tw

CAMECA Instruments, Inc. USA • cameca.us-sales@ametek.com

 **CAMECA**[®]
SCIENCE & METROLOGY SOLUTIONS

 **AMETEK**[®]
MATERIALS ANALYSIS DIVISION

www.cameca.com